

# LECO – Your Source for Total Analytical Solutions




### TGA801 Thermogravimetric Analyser

- Determines moisture and ash or other gravimetric methods in various organic samples
- Expanded temperature control (up to 1,000 °C) with variable atmospheres and ramp rates
- Automated batch thermogravimetric analysis of up to 19 samples



### AC Series Isoperibol Calorimeters

- Accurate calorific measurements
- Ergonomic water-measuring/combustion vessel-filling station
- Analysis time in as little as five minutes
- Meets or exceeds ASTM requirements
- 6,000 to 15,000 BTU/lb. for 1 gram sample



### Organic Consumables

Get the best results from your LECO instrument by using genuine LECO consumables. Visit [eu.leco.com](http://eu.leco.com) for featured items, specials, and ordering information (Form number 203-828).

## 832 Series Sulfur/Carbon Analysis by Combustion



ISO-9001:2015 Q-994 | LECO, Cornerstone are a registered trademark of LECO Corporation.

**LECO Corporation**  
3000 Lakeview Avenue | St. Joseph, MI 49085  
Phone: 269-985-5496  
[info@leco.com](mailto:info@leco.com) | [www.leco.com](http://www.leco.com)

**LECO Europe**  
[eu.leco.com](http://eu.leco.com)



# 832 Series: Sulfur / Carbon Analysis by Combustion

LECO's 832 Series redefines the way you determine carbon and sulfur in coal, coke, cement, soil, fuel oils, and other organic materials. Using extensive customer feedback and innovative engineering, our latest instrumentation takes advantage of an on-board software platform powered by a touch-screen interface, increasing usability and operational control without using valuable bench space.


Additional features – such as a high-efficiency combustion furnace, improved IR cell design, rugged design, and a robust, rugged design – make the 832 a valuable resource for any laboratory needing fast and accurate analysis of sulfur and/or carbon in organic materials.



# User-friendly CORNERSTONE® Brand Software

LECO's exclusive Cornerstone brand software with touch-screen interface enables the user to have complete access to analysis control, method settings, diagnostics, reporting, and more in a highly organised and immersive environment. Designed through a collaboration of customer feedback and innovative engineering,

Cornerstone allows the user to conduct all of their day-to-day operations within a single analysis screen designed for speed and ease of use. Our innovative grouping of sample data into sets and replicates simplifies the data output and automatically calculates relevant statistics, alleviating the need for additional data processing.




Analysis



# Software Features and Benefits

- The software is divided into four main sections – Analysis, Diagnostics, Settings, and Instrument – for simplified navigation and organisation
- Toolbars, sliders, and drop-down menus make it easy to set parameters for calibration and data processing
- The software also includes real-time monitoring of ambient parameters, with fully animated system diagrams
- Advanced interactive diagnostic features include a thorough digital on-board manual, maintenance animations, photo illustrations, and screen captures that quickly provide the direction needed without having to refer to multiple manuals
- CORNERSTONE® also supports a multilingual interface, user permissions, extended data archiving and filtering, compatibility with various Laboratory Information Management Systems (LIMS), and flexible reporting capabilities





Diagnostics > Ambients



Instrument > Furnace




Settings > Calibration




An optional Cornerstone Mobile application feature enables remote viewing of the instrument software from a smart phone, tablet, or PC. It can also be programmed to set automatic notifications from the instrument against predefined software conditions using e-mail, text message, or the Cornerstone Mobile application.

# Instrument Highlights and Features



## High-Efficiency Furnace System Design

- Lower electrical requirement resulting in lower operating costs
- Intelligent furnace control optimises furnace reliability and extends the heating element and ceramics lifetime
- Horizontal high temperature (up to 1,550 °C with optional HT model) concentric ceramic furnace tube with ceramic lance directing oxygen directly onto the sample ensuring rapid analysis times, complete combustion, and recovery
- Improved gas flow control reduces ash debris restrictions, increases furnace lifetime, and instrument reliability
- Modified combustion assembly provides furnace isolation from the atmosphere, resulting in low blank value and improved precision



## Improved Solid-State IR Cell Design

- Temperature stabilised construction provides increased isolation from ambient temperature fluctuations and other environmental susceptibilities
- Optimised emitter control and detection circuitry improves the IR cell lifetime and long-term stability, resulting in superior accuracy and precision
- Individual wide-range IR detection for both carbon and sulfur simplifies set-up and calibration for the standard S/C/SC832 models
- Optional dual-range (DR) detection for sulfur within S/SC832DR models have the widest sulfur range capability for the most demanding applications



Improved Design and Accessibility for Maintenance Areas

- To enhance safety and convenience, the dual reagent tubes and rotameters are located on the front side of the instrument behind the hinged front cabinet door
- Open internal cabinet access to pump with quick release features enables fast and simple preventive maintenance routines, ensuring a robust and reliable instrument with superior uptime
- Top-mounted heating element assemblies result in simple and fast access, and allow for single element replacement when needed, lowering maintenance costs



Additional Features and Benefits

- Optional Autoloader (L) with 100-sample capacity simplifies operator workflow and increases productivity
- Optional High Temperature (HT) model supports furnace temperatures up to 1,550 °C, eliminating the need for combustion accelerators for refractory sample analysis
- Optional boom-mounted touch-screen interface package (M) promotes an ergonomic workspace while reducing system bench space requirements
- Optional Performance Package (P) adds an instrument leak check feature, which provides a valuable diagnostic tool and increases the reliability of the system; and an electronic back pressure control, which provides isolation to changes in atmospheric pressure and results in improved long-term stability
- Dual range (DR) S/SC832DR models contain the performance package features

Models

S832	Sulfur Determinator
C832	Carbon Determinator
SC832	Sulfur and Carbon Determinator
S832DR	Dual Range Sulfur Determinator
SC832DR	Dual Range Sulfur with Carbon Determinator
S832HT	High Temperature Dual Range Sulfur Determinator
SC832HT	High Temperature Dual Range Sulfur with Carbon Determinator

Options

Monitor (M)	Instrument mounted touch-screen and six-axis boom mount
Performance Pack (P)	Adds an instrument leak check and electronic back pressure control (All DR instrument models contain performance pack)
Autoloader (L)	100-capacity sample Autoloader

NOTE: Multiple configurations of options are available.  
Please contact your local LECO Sales Engineer for more details.